




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz reaction vessel, cylindrical, with laboratory flange LF 200, 10,000 ml capacity, outer/inner Ø 215/201 mm, height approx. 450 mm, with round rim**

**€1,104.05  
plus VAT &  
Shipping**



## Description

---

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube.

## Additional Information

---

No.	LZ-1662
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3665
EAN	4058072314002
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	200 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	10 l
Form	cylindrical neck
Height	450 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	200
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

