



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction vessel, cylindrical, with laboratory flange LF 200, 20,000 ml capacity, outer/inner Ø 315/301 mm, height approx. 465 mm, with round bottom

**€1,597.02
plus VAT &
Shipping**

Product Images



Description

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube.

Additional Information

No.	LZ-1663
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3666
EAN	4058072314019
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	300 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	20 l
Form	cylindrical neck
Height	465 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	200
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

