



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Preparation tubes 5 ml, grad., conical, with lid, 250 pcs./pack

€79.00
plus VAT &
Shipping

Product Images



Description

Conical preparation tubes made of PP. Ideal for preparing light-sensitive biological or chemical solutions (black tubes), also for microbiological incubation and pelleting. Tubes are graduated in 1 ml sections and can be labeled on frosted areas on the lid and side.

Temperature resistant from -80°C to 121°C.

Centrifugable up to 5000 x g.

Dimensions: 18 x 60 mm (Ø x H)

Additional Information

No.	7-4730
Manufacturer (Brand)	neoLab
EAN	4058072041144
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	18 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	5 ml
sterile	No
Height	60 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	121 °C
TBST MIN	-80 °C
Reaction vessels and test tubes Properties	centrifugable surface-optimized with closure/lid with labeling field graduated
Floor outside	conical
Floor inside	conical
Closure / lid properties	closed (attached)
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Type Tube/ Micro Vessel	Reaction vessels Preparation tubes
centrifugable up to	5000 g

