



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® reaction tubes 1.5 ml, with tamper-evident closure 500 pcs/pack, green

**€27.50
plus VAT &
Shipping**

Product Images



Description

Reaction vessels with tamper-evident closure for sensitive samples, e.g. in forensic medicine. During sample preparation, it is ensured that the contents of the reaction vessel have already been processed when the lid is closed and the tamper-evident closure is intact.

Can be used at -180°C to $+121^{\circ}\text{C}$. Centrifugable up to 20,000 RZB.

Material PP with very low adhesion, labeling field.

Graduated at 0.1, 0.5, 1 and 1.5 ml.

Additional Information

No.	7-4134
Manufacturer (Brand)	neoLab
EAN	4058072039011
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	green
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
Type closure/ lid	Push-on/plug-in closure original closure
TBST MAX	121 °C
TBST MIN	-85 °C
Reaction vessels and test tubes Properties	centrifugable surface-optimized with closure/lid with labeling field graduated
Floor outside	conical
Floor inside	conical
Type Tube/ Micro Vessel	Reaction vessels
Closure / lid properties	closed (attached)
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
centrifugable up to	15000 g

