



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® reaction tubes w. lid 0.5 ml, violet, 1000
pcs./pack**

**€50.20
plus VAT &
Shipping**

Product Images



Description

Colored reaction vessels made of PP. Tightly closing lid, tight fit, easy to open. Lid membrane is uniform thickness and large to total Ø. The different colors guarantee easy differentiation.

Additional Information

No.	1-6165
Manufacturer (Brand)	neoLab
EAN	4058072036980
autoclavable at 121°C	Yes
diameter	7 mm
Color	purple
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	No
Height	30 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	121 °C
TBST MIN	-80 °C
Reaction vessels and test tubes Properties	centrifugable surface-optimized with closure/lid with labeling field graduated
Floor outside	conical
Floor inside	conical
Type Tube/ Micro Vessel	Reaction vessels
Closure / lid properties	closed (attached)
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
centrifugable up to	20000 g

