



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Reaction vessels w. lid 0.5 ml, orange,  
1000 pcs./pack**

**€50.20  
plus VAT &  
Shipping**

## Product Images



## Description

---

Orange reaction tubes made of PP. Tightly closing lid, tight fit, easy to open. Lid membrane is uniform thickness and large to total Ø. The different colors guarantee easy distinction.

## Additional Information

---

No.	1-6166
Manufacturer (Brand)	neoLab
EAN	4058072036997
autoclavable at 121°C	Yes
diameter	7 mm
Color	orange
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	No
Height	30 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	121 °C
TBST MIN	-80 °C
Reaction vessels and test tubes Properties	centrifugable surface-optimized with closure/lid with labeling field graduated
Floor outside	conical
Floor inside	conical
Type Tube/ Micro Vessel	Reaction vessels
Closure / lid properties	closed (attached)
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
centrifugable up to	20000 g

