



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz reaction vessel, spherical, with  
temperature control jacket, with laboratory  
flange LF 100, 4,000 ml capacity, outer/inner Ø  
236/206 mm, height approx. 285 mm, spigot GL  
18**

**€1,461.08  
plus VAT &  
Shipping**



## Description

---

With laboratory flange LF, without pitch, with groove. 2 nozzles with GL thread. Made of DURAN® tube.

## Additional Information

---

|   |                           |
|---|---------------------------|
| No.   | LZ-1747                   |
| Manufacturer (Brand)                                    | Lenz                      |
| old neoLab article no.                                  | B-3750                    |
| EAN   | 4058072314859             |
| Transport temperature                                   | Room temperature          |
| autoclavable at 121°C                                   | Yes                       |
| diameter  | 236 mm                    |
| Color   | transparent (clear)       |
| Material  | Borosilicate glass        |
| Volume  | 4 l                       |
| Form  | spherical                 |
| Height  | 295 mm                    |
| Connection General / Tempering jacket / Cooler / Vacuum | GL 18                     |
| Laboratory plan grinding                                | 100                       |
| Type glass jar  | Reaction vessel spherical |
| Type ground   | Laboratory plan grinding  |

