



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction vessel, spherical, with temperature control jacket and drain valve, with laboratory flange LF 100, 4,000 ml capacity, outer/inner Ø 236/206 mm, height approx. 395 mm, with 2 nozzles GL 18, valve 10 mm

**€2,054.48
plus VAT &
Shipping**



Description

With laboratory flange LF, without pitch, with groove. 2 spigots (with thread GL) Made of DURAN® pipe

Additional Information

No.	LZ-1798
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3801
EAN	4058072315368
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	236 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	4 l
Form	spherical
Height	395 mm
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Laboratory plan grinding	100
Type glass jar	Reaction vessel spherical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

