



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction vessel, spherical, with temperature control jacket and drain valve, with laboratory flange LF 60, 500 ml capacity, outer/inner Ø 131/105 mm, height approx. 280 mm, with 2 nozzles KF 15, valve 10 mm

**€1,818.82
plus VAT &
Shipping**



Description

With laboratory flange LF, without pitch, with groove. 2 spigots (with ball flanges KF) Made of DURAN®pipe

Additional Information

No.	LZ-1802
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3805
EAN	4058072315405
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	131 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	0.5 l
Form	spherical
Height	280 mm
Connection General / Tempering jacket / Cooler / Vacuum	Ball flange 15
Laboratory plan grinding	60
Type glass jar	Reaction vessel spherical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

