



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction vessel, spherical, with temperature control jacket and drain valve, with laboratory flange LF 150, 6,000 ml capacity, outer/inner Ø 280/236 mm, height approx. 455 mm, with 2 nozzles KF 15, valve 10 mm

**€2,513.30
plus VAT &
Shipping**



Description

With laboratory flange LF, without pitch, with groove. 2 spigots (with ball flanges KF) Made of DURAN®pipe

Additional Information

No.	LZ-1807
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3810
EAN	4058072315450
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	280 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	6 l
Form	spherical
Height	455 mm
Connection General / Tempering jacket / Cooler / Vacuum	Ball flange 15
Laboratory plan grinding	150
Type glass jar	Reaction vessel spherical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

