




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® reaction tubes black, 1.5 ml, 500
pcs./pack**

**€32.90
plus VAT &
Shipping**

Product Images



Description

Black reaction tubes made of PP

For storage of light-sensitive samples and for other applications with light-protected storage.

The 1.5 ml tube is graduated at 0.1, 0.5, 1.0 and 1.5 ml.

Tapered, temperature resistant from -80°C to +121°C.

Centrifugable up to 15,000 x g.

Material: PP, autoclavable

Additional Information

No.	7-7590
Manufacturer (Brand)	neoLab
EAN	4058072041670
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	Black
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	121 °C
TBST MIN	-80 °C
Reaction vessels and test tubes Properties	centrifugable surface-optimized with closure/lid with labeling field graduated for light-sensitive samples
Floor outside	conical
Floor inside	conical
Closure / lid properties	closed (attached)
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	15000 g

