



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® reaction tubes 0.5 ml, conical, with
screw cap, 500 pcs./pack**

**€138.20
plus VAT &
Shipping**

Product Images



Description

Reaction vessels made of polypropylene, conical.

With screw cap and EPR sealing ring (=ethylene propylene rubber) for secure, leak-free closure.

Handle rim slightly ribbed for easier opening, especially when gloves are worn.

Additional Information

No.	3-1942
Manufacturer (Brand)	neoLab
EAN	4058072018542
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	7 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	No
Height	30 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable surface-optimized with closure/lid with seal can be used as cryotube
Floor outside	conical
Floor inside	conical
Material gasket / septum	Ethylene propylene rubber (EPR)
Type Tube/ Micro Vessel	Cryotubes Reaction vessels
Closure / lid properties	closed
suitable for cooling in	Nitrogen tank
centrifugable up to	20000 g

