



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® reaction tubes graduated, conical, 1.5 ml, 1000 pcs./pack

Special Price

~~€99.00~~ was

**€148.50**

**plus VAT &  
Shipping**



## Description

---

Reaction vessels, graduated, with white labeling field.  
Self-standing version with stand ring.

Please order caps separately.

## Additional Information

---

No.	7-4527
Manufacturer (Brand)	neoLab
EAN	4058072040420
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized with labeling field graduated can be used as cryotubes
Floor outside	conical
Floor inside	conical
suitable for cooling in	Nitrogen tank
Type Tube/ Micro Vessel	Reaction vessels Cryogenic tubes
centrifugable up to	17000 g

