



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® reaction tubes graduated, self-
standing, 2.0 ml, 1000 pcs./pack**

**€169.90
plus VAT &
Shipping**

Product Images



Description

Reaction vessels, graduated, with white labeling field.
Self-standing version with stand ring.

Please order caps separately.

Additional Information

No.	7-4528
Manufacturer (Brand)	neoLab
EAN	4058072040437
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	2 ml
sterile	No
Height	42 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized with labeling field graduated can be used as cryotubes
Floor outside	conical. pointed . with stand ring
Floor inside	conical
suitable for cooling in	Nitrogen tank
Type Tube/ Micro Vessel	Reaction vessels Cryogenic tubes
centrifugable up to	17000 g

