



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® reaction tubes ungraduated, 1.5 ml,
self-standing, 1000 pcs./pack**

**€98.50
plus VAT &
Shipping**

Product Images



Description

Reaction vessels, ungraduated

- Self-standing version with stand ring

In resealable bag.

Please order closures separately.

Additional Information

No.	7-4521
Manufacturer (Brand)	neoLab
EAN	4058072040369
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized can be used as cryotubes
Floor outside	conical. with stand ring
Floor inside	conical
TBST MAX	121 °C
TBST MIN	-196 °C
Type Tube/ Micro Vessel	Reaction vessels Cryogenic tubes
suitable for cooling in	Nitrogen tank
centrifugable up to	17000 g

