



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® reaction tubes ungraduated, 2.0 ml,  
conical, 1000 pcs./pack**

**€98.50  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

### Reaction vessels, ungraduated

- conical design

In resealable bag.

Please order closures separately.

## Additional Information

---

No.	7-4524
Manufacturer (Brand)	neoLab
EAN	4058072040390
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	2 ml
sterile	No
Height	42 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized can be used as cryotubes
Floor outside	conical
Floor inside	conical
Type Tube/ Micro Vessel	Reaction vessels Cryogenic tubes
suitable for cooling in	Nitrogen tank
centrifugable up to	17000 g

