



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand reaction vessels 1.5 ml, transparent w. lid,
500 pcs./pack**

€22.25 As low as
€21.60
**plus VAT &
Shipping**



Description

Reaction vessels, PP, 1.5ml with lid

These reaction tubes made of PP with a volume of 1.5 ml incl. lid are free of PCR inhibitors, perfectly sealing to protect against evaporation and are characterized by high transparency and a pierceable lid membrane.

German quality goods - centrifugable up to 20000 RZB, frosted labeling field, with graduation.

Additional Information

No.	BN-2138
Manufacturer (Brand)	Brand
old neoLab article no.	E-1028
EAN	4058072213626
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	121 °C
TBST MIN	-80 °C
Reaction vessels and test tubes Properties	centrifugable surface-optimized with closure/lid with labeling field graduated
Floor outside	conical
Floor inside	conical
Closure / lid properties	closed (attached)
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	15000 g

