



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® reaction tubes 0.5 ml, transparent w. lid, 1000 pcs./pack

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Reaction vessels, PP

transparent, German quality goods.

Additional Information

No.	E-1030
Manufacturer (Brand)	neoLab
EAN	4058072062200
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	7 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	No
Height	30 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	121 °C
TBST MIN	-80 °C
Reaction vessels and test tubes Properties	centrifugable surface-optimized with closure/lid with labeling field graduated
Floor outside	conical
Floor inside	conical
Closure / lid properties	closed (attached)
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	15000 g

