



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® Reaction vessels w. screw thread, 1.5 ml, conical, 500 pcs./pack**

**€44.30**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

**Reaction vessels (PP) with screw thread, transparent**, with uniform wall thickness.

Packing: 500 pieces in a bag.

## Additional Information

---

No.	1-6186
Manufacturer (Brand)	neoLab
EAN	4058072037086
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized can be used as cryotubes
Floor outside	conical
Floor inside	conical
suitable for cooling in	Nitrogen tank
Type Tube/ Micro Vessel	Reaction vessels Cryogenic tubes
centrifugable up to	20000 g

