



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Reaction vessels w. screw thread, 2.0 ml, self-standing, 500 pcs./pack

€44.30
plus VAT &
Shipping

Product Images



Description

Reaction vessels (PP) with screw thread, transparent, with uniform wall thickness.

Packaging: 500 pieces in a bag.

Additional Information

No.	1-6190
Manufacturer (Brand)	neoLab
EAN	4058072042813
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	2 ml
sterile	No
Height	42 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized can be used as cryotubes
Floor outside	conical. with stand ring
Floor inside	conical
Type Tube/ Micro Vessel	Reaction vessels Cryogenic tubes
suitable for cooling in	Nitrogen tank
centrifugable up to	20000 g

