



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction vessel, cylindrical, with temperature control and vacuum jacket and drain valve, with laboratory flange LF 100, 3,000 ml, outer/inner Ø 250/150 mm, height approx. 310 mm, with 2 nozzles KF 25

**€2,426.65
plus VAT &
Shipping**

Product Images



Description

With laboratory flange LF, with graduation and groove. 2 spigots; ball flanges. Made of DURAN® pipe.

Additional Information

No.	LZ-1819
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3822
EAN	4058072315573
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	250 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	3 l
Form	cylindrical
Height	310 mm
Reaction vessels and test tubes Properties	graduated
Connection General / Tempering jacket / Cooler / Vacuum	Ball flange 25
Color graduation	black
Laboratory plan grinding	100
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding
Type valve	Throttle valve with locking screw

