



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Safety-Cap Reaction Tubes transparent, sterile, 0.5 ml 50 pieces/pack

€42.60
plus VAT &
Shipping

Product Images



Description

The **sterile Safety-Cap reaction tubes** feature a stable lid lock that prevents uncontrolled opening of the tubes, e.g. during incubation. Especially when working with valuable samples as well as radioactive substances or DNA, a secure closure is essential. The reaction vessels are IVD certified and therefore also suitable for biological samples in in-vitro diagnostics.

- autoclavable in opened condition at 121 °C, 20 min
- 0.5 ml to 30,000 g centrifugable
- DNase-/RNase-, PCR inhibitors, endotoxin-, DNA- and ATP-free
- high chemical resistance

Additional Information

No.	E-0180
Manufacturer (Brand)	neoLab
EAN	4058072160951
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	8 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	Yes
Height	31 mm
Type closure/ lid	Pressure cap
TBST MAX	121 °C
TBST MIN	0 °C
Reaction vessels and test tubes Properties	centrifugable with cap/lid RNase-free DNase-free DNA-free
Floor outside	conical
Floor inside	conical
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator freezer/ freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	30000 g

