



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Safety-Cap Reaction Vials transparent 1.5 ml ( Pack 1000 pcs )

**€43.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The **Safety-Cap reaction tubes** feature a stable lid lock that prevents uncontrolled opening of the tubes, e.g. during incubation. Especially when working with valuable samples as well as radioactive substances or DNA, a secure closure is essential.

- autoclavable when open at 121 °C, 20 min
- 1.5 ml to 30,000 g centrifugable
- DNase-/RNase-, PCR inhibitors, endotoxin-, DNA- and ATP-free
- high chemical resistance

## Additional Information

No.	E-1015
Manufacturer (Brand)	neoLab
EAN	4058072112530
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	40 mm
Type closure/ lid	Pressure cap
TBST MAX	121 °C
TBST MIN	0 °C
Reaction vessels and test tubes Properties	centrifugable with cap/lid RNase-free DNase-free DNA-free
Floor outside	conical
Floor inside	conical
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator freezer/ freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	30000 g

