



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® neoScrew Microtubes, low wetting,  
self-standing, 1.5 ml, 1000 pcs/pack**

**€110.30  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Reaction vessels made of natural-colored special polypropylene for low-wetting surfaces, ungraded. Ideal for nucleic acid or protein work.

Please order screw caps separately.

## Additional Information

---

No.	7-4713
Manufacturer (Brand)	neoLab
EAN	4058072041045
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized can be used as cryotubes
Floor outside	conical. with stand ring
Floor inside	conical
Type Tube/ Micro Vessel	Reaction vessels Cryogenic tubes
suitable for cooling in	Nitrogen tank
centrifugable up to	17000 g

