



neolab Migge GmbH

Rischerstr. 7-9

69123 Heidelberg

Deutschland

+49 (0)6221 / 8442-44

<https://www.neolab.de>

Umsatzsteuer-Identifikationsnummer:

DE 143 450 657



neolab® neoScrew microtubes, low wetting, conical bottom, 1.5 ml, 1000 pcs/pack

€110.30

plus VAT & Shipping



Description

Reaction vessels made of natural-colored special polypropylene for low-wetting surfaces, ungraded.

Ideal for nucleic acid or protein work.

Please order screw caps separately.

Additional Information

No.	7-4714
Manufacturer (Brand)	neoLab
EAN	4058072041052
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized can be used as cryotubes
Floor outside	conical
Floor inside	conical
suitable for cooling in	Nitrogen tank
Type Tube/ Micro Vessel	Reaction vessels Cryogenic tubes
centrifugable up to	17000 g

