



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® neoScrew Micro Tubes, brown, self-
standing, 0.5 ml, 1000 pcs./pack**

**€116.10
plus VAT &
Shipping**

Product Images



Description

Reaction vessels made of dark brown PP, ideal for light-sensitive samples.
Not graduated. With screw thread.

Please order screw caps separately.

Additional Information

No.	7-4590
Manufacturer (Brand)	neolab
EAN	4058072040826
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	6 mm
Color	brown
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	No
Height	31 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized can be used as cryotubes
Floor outside	conical. with stand ring
Floor inside	conical
suitable for cooling in	Nitrogen tank
Type Tube/ Micro Vessel	Reaction vessels Cryogenic tubes
centrifugable up to	17000 g

