



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® neoScrew Microtubes 2.0 ml, sterile
with cap, self-standing, 10x50 pcs/pack**

**€177.90
plus VAT &
Shipping**

Product Images



Description

Sterile reaction vessels made of PP with thread and cap,
tightly closing due to **high-quality O-ring seal** .
Ideal for cryogenic work. Ungraded.

Sterile packed in 10 bags of 50 pieces.

Additional Information

No.	7-4578
Manufacturer (Brand)	neoLab
EAN	4058072040758
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	2 ml
sterile	Yes
Height	47 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized with closure/lid with seal can be used as cryotube
Floor outside	flat conical. with stand ring
Floor inside	conical. pointed
Material gasket / septum	Ethylene propylene rubber (EPR)
Type Tube/ Micro Vessel	Cryotubes Reaction vessels
Closure / lid properties	closed
suitable for cooling in	Nitrogen tank
centrifugable up to	17000 g

