



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® neoScrew-Microtubes 0.5 ml self-
standing, sterile, cap+flap, 10x50 pcs./pack**

**€227.80
plus VAT &
Shipping**

Product Images



Description

Sterile reaction vessels (PP) with cap with tab, tightly closing due to high-quality O-ring seal.

Ideal for cryo work.

Ungraded, sterile packed, in 10 bags of 50 pieces.

Additional Information

No.	7-4700
Manufacturer (Brand)	neoLab
EAN	4058072040932
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	Yes
Height	47 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized with closure/lid with seal can be used as cryotube
Floor outside	flat conical. with stand ring
Floor inside	conical
Material gasket / septum	Ethylene propylene rubber (EPR)
Closure / lid properties	closed (attached)
suitable for cooling in	Nitrogen tank
Type Tube/ Micro Vessel	Cryotubes Reaction vessels
centrifugable up to	17000 g

