



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® neoScrew-Microtubes 1,5 ml conical,  
sterile grad., cap+flap, 10x50 pcs./pack**

**€256.20  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Sterile reaction vessels made of PP with cap with tab.

Tightly closing due to O-ring seal.

Ideal for cryogenic work.

Graduated, with white labeling field.

Sterile packed in 10 bags of 50 pieces.

## Additional Information

No.	7-4708
Manufacturer (Brand)	neoLab
EAN	4058072041007
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	Yes
Height	47 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized with closure/lid with seal with labeling field graduated can be used as cryotube
Floor outside	conical
Floor inside	conical
Material gasket / septum	Ethylene propylene rubber (EPR)
Type Tube/ Micro Vessel	Cryotubes Reaction vessels
Closure / lid properties	closed (attached)
suitable for cooling in	Nitrogen tank
centrifugable up to	17000 g

