



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Beckman centrifuge bottles, 250 ml, 62
x 120 mm, PC, 6 pcs/pack**

**€191.75
plus VAT &
Shipping**

Product Images



Description

Beckman-Coulter High-Performance Centrifuge Bottles for highest loads up to 110000 x g.
All bottles complete with tight-fitting caps.

Additional Information

No.	C-1132
Manufacturer (Brand)	neoLab
EAN	4058072052980
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Speed MAX	5200 rpm
diameter	62 mm
Color	transparent (crystal clear)
Volume	250 ml
sterile	No
Height	120 mm
Type closure/ lid	Screw cap
TBST MAX	125 °C
TBST MIN	0 °C
Floor outside	flat
Type Tube/ Micro Vessel	Centrifuge bottles
centrifugable up to	38400 g

