



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® high-speed centrifuge tube 30 ml

**€71.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

High-speed centrifuge glass made of reinforced borosilicate glass.  
Virtually unbreakable due to special chemical treatment.  
Can withstand up to 13,100 RZB when used correctly.  
Round-bottomed, without rim.  
Maximum temperature: 300°C  
Autoclavable at 121°C (15 min)

## Additional Information

---

No.	2-9211
Manufacturer (Brand)	neoLab
EAN	4058072010461
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	24 mm
Color	transparent (clear)
Volume	30 ml
sterile	No
Height	106 mm
Type closure/ lid	without closure
TBST MAX	300 °C
TBST MIN	-70 °C
Floor outside	around
Material tube	Borosilicate glass
Type Tube/ Micro Vessel	Centrifuge tubes
Wall thickness	normal
centrifugable up to	12000 g

