



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Falcon® tubes conical, PS, in styrofoam rack, 15 ml, 10x50 pcs.

**€202.44
plus VAT &
Shipping**

Product Images



Description

Centrifuge tubes conical made of polystyrene for highest optical requirements.

- max. centrifugation 1,800 g
- tight-fitting, self-sealing screw cap
- white graduation from 1.5 ml to 14.5 ml in 0.5 ml steps
- white labeling field
- sterile, gamma irradiated

PU: 500 pieces (10 x 50)

Order No. C-3081

In polystyrene rack

Additional Information

No.	GF-0056
Manufacturer (Brand)	Falcon
old neoLab article no.	C-3081
EAN	4058072256173
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	17 mm
Color	transparent (crystal clear)
Volume	15 ml
sterile	Yes
Height	120 mm
Type closure/ lid	Screw cap
TBST MAX	60 °C
TBST MIN	4 °C
Floor outside	conical
Color graduation	white
Type of sterility	radiation sterile
Division	0.5 ml
Format/ Packing	in bag
Material tube	Polystyrene (PS)
Material closure/ lid	Polystyrene (PS)
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	1800 g

