



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Moonlab® Centrifuge tubes 50 ml, conical,
clear, flat cap, sterile, PP, 20 x 25 pcs, in bag**

Special Price

€84.00 was

€120.00

**plus VAT &
Shipping**

Product Images



Description

50 ml conical centrifuge tubes made of PP with screw cap made of PE.

- centrifugable up to 8.600 rpm for 15 ml and up to 9.000 rpm for 50 ml
- Temperature resistance: -20 °C to 121 °C
- pyrogen-free
- sterile (gamma radiation)
- graduated with labeling field
- ideal for all standard applications in the life science sector

The tubes are available in racks or bags.

PU: 500 pcs./pack

Order-No. 4-0003

in bag, underpacked 20 x 25 pcs

Additional Information

No.	4-0003
Manufacturer (Brand)	Moonlab
EAN	4058072084905
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	30 mm
Color	transparent (clear)
Volume	50 ml
sterile	Yes
Height	121 mm
Type closure/ lid	Screw cap
Type of sterility	radiation sterile
Floor outside	conical
Color graduation	white
Format/ Packing	in bag
Material tube	Polypropylene (PP)
Material closure/ lid	High-density polyethylene (HDPE)
TBST MAX	121 °C
TBST MIN	-20 °C
Division	1 ml
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	9400 g

