



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene Oak Ridge Centrifuge Tubes for High Speed Applications 50mL Case of 10

€168.24 As low
as
€135.40
plus VAT &
Shipping

Product Images



Description

Nalgene™ Oak Ridge Centrifuge Tubes for High-Speed Applications, 50 mL Thermo Scientific Nalgene Oak Ridge centrifuge tubes made of FEP for high-speed applications are excellent for high-speed centrifugation of aggressive chemicals. Closure Material: ETFE Capacity (metric): 50 mL Number per box: 10 Number per pack: 2 Autoclavable: Yes Graduations: No Shape: Round Type: centrifuge tubes Max. Force: 50,000 × g Closure Type: Screw Cap Writing Area: No Unit Size: Case of 10

Additional Information

No.	CD-0752
Manufacturer (Brand)	Nalgene
VGKL number	238731143
EAN	4058072171544
Gross weight	3.24
Transport temperature	Room temperature
Material	Polyfluoroethylene propylene (FEP)
Volume	50 ml
Country of origin	USA
Packaging volume	0.00200000 ccm
Customs tariff number	39269097

