



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene PPCO centrifuge tubes with round bottom for high-speed applications

€69.00 As low
as
€55.60
plus VAT &
Shipping

Product Images

Description

Nalgene™ PPCO round-bottom centrifuge tubes for high-speed applications Thermo Scientific™ Nalgene™ PPCO round-bottom centrifuge tubes for high-speed applications have uniform wall thickness and diameter, fit better in rotors, and perform better. Capacity (metric): 50 mL Quantity per case: 100 Quantity per pack: 10 Autoclavable: yes Graduations: No Shape: Round Type: Centrifuge Tube Max. Force: 50,000 x g Writing Area: No Unit Size: Case of 100

Additional Information

No.	CD-0748
Manufacturer (Brand)	Nalgene
EAN	4058072174460
Gross weight	2.16
Transport temperature	Room temperature
Material	Polypropylene (PP)
Volume	50 ml
Country of origin	USA
Packaging volume	0.03600000 cc
Customs tariff number	39269097

