



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Falcon® round bottom tubes, PS, 14 ml, 2-pos.  
closure, sterile, 20x25 pcs.**

**€126.71  
plus VAT &  
Shipping**

## Product Images



## Description

---

Centrifuge tubes with round bottom made of polystyrene for the use of various test preparations.

- max. centrifugation 1.400 g
- highest transparency
- tightly closing 2-position plug-in cap for pressure or gas exchange with simultaneous lid lock
- white graduation from 1.5 ml to 14.5 ml in 0.5 ml steps
- Graduation printed or embossed
- sterile, gamma irradiated

**PU:** 500 pieces (20 x 25)

## Additional Information

No.	GF-0044
Manufacturer (Brand)	Falcon
old neoLab article no.	C-3085
EAN	4058072256111
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	17 mm
Color	transparent (crystal clear)
Volume	14 ml
sterile	Yes
Height	100 mm
Type closure/ lid	Push-on/plug-in closure with pressure/gas exchange
TBST MAX	60 °C
TBST MIN	4 °C
Floor outside	around
Color graduation	white
Type of sterility	radiation sterile
Division	0.5 ml
Format/ Packing	in bag
Material tube	Polystyrene (PS)
Material closure/ lid	Polystyrene (PS)
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	1400 g

