



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Falcon® round bottom tubes, PS, 5 ml, 2-pos.
closure, sterile, 20x25 pcs.**

**€120.40
plus VAT &
Shipping**

Product Images



Description

Centrifuge tubes with round bottom made of polystyrene for the use of various test preparations.

- max. centrifugation 1.400 g
- highest transparency
- tightly closing 2-position plug cap for pressure or gas exchange with simultaneous lid lock
- white graduation from 1.5 ml to 14.5 ml in 0.5 ml steps
- Graduation printed or embossed
- sterile, gamma irradiated

PU: 500 pieces (20 x 25)

Additional Information

No.	GF-0045
Manufacturer (Brand)	Falcon
old neoLab article no.	C-3082
EAN	4058072256128
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	12 mm
Color	transparent (crystal clear)
Volume	5 ml
sterile	Yes
Height	75 mm
Type closure/ lid	Push-on/plug-in closure with pressure/gas exchange
Type of sterility	radiation sterile
Floor outside	around
Color graduation	white
Format/ Packing	in bag
Material tube	Polystyrene (PS)
Material closure/ lid	Polystyrene (PS)
TBST MAX	60 °C
TBST MIN	4 °C
Division	0.5 ml
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	1400 g

