



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® urine tubes made of PS with lid 500
pcs/pack**

**€110.50
plus VAT &
Shipping**

Product Images



Description

neoLab urine tubes made of polystyrene with matching tight-fitting lids

Urine tubes are made of pure polystyrene and are free of chemicals, metals or additives that could contaminate specimens.

Thick walled; graduated 1/4, 1/2, 1, 2 and every 2 ml up to 12 ml.

Tubes can be safely centrifuged at speeds up to 2000 g.

They are flared at the top to facilitate filling and drip-free pouring. Designed to allow the use of miniature urinometers and reagent test strips that require only 1/4 or 1/2 ml of sample.

Additional Information

No.	2-1023
Manufacturer (Brand)	neoLab
EAN	4058072155612
autoclavable at 121°C	No
diameter	21 mm
Color	transparent (clear)
Volume	12 ml
sterile	No
Height	105 mm
Type closure/ lid	Screw cap
Floor outside	conical
Format/ Packing	in bag
Material tube	Polystyrene (PS)
Material closure/ lid	Polypropylene (PP)
Type Tube/ Micro Vessel	Centrifuge tube
centrifugable up to	2000 g

