



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hecht centrifuge tubes round, Duran, 100 x 16 mm, 10 pcs/pack

€6.57 As low as

€4.33

**plus VAT &
Shipping**

Product Images

Description

Without rim, made of borosilicate glass 3.3, German quality goods.

Additional Information

No.	WQ-0625
Manufacturer (Brand)	Hecht
old neolab article no.	E-1078
EAN	4058072443481
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	16 mm
Color	transparent (clear)
Volume	12 ml
sterile	No
Height	100 mm
TBST MAX	300 °C
TBST MIN	-70 °C
Floor outside	around
Format/ Packing	individually packed in carton/ dispenser
Material tube	Borosilicate glass
Type Tube/ Micro Vessel	Centrifuge tubes
Wall thickness	strong-walled

