



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific EZFlip centrifuge tubes conical, 15 ml, sterile, 500 pcs/pack

€282.00
plus VAT &
Shipping

Product Images



Description

Thermo Scientific™ Nunc™ EZFlip™ conical centrifuge tubes can be opened with one hand thanks to the proprietary system with attached closures open and close with one hand. The tubes have an ergonomic design and are compatible with all fixed-angle or swing-out rotors. The tubes are sterilized and supplied in bulk (in a pouch).

- The attached, leak-proof cap prevents cross-contamination and loss of the cap.
- Conical bottom, sterile
- Clarified polypropylene is chemical resistant and withstands high centrifugal forces
- 15 ml withstands accelerations up to 8,500 x g
- White title block and flat lid for sample identification
- Printed volume markings and graduations
- Standard dimensions

Additional Information

No.	FE-0674
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	7-5008
EAN	4058072254094
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Volume	15 ml
sterile	Yes
Type closure/ lid	Flip closure (one-handed operation)
Type of sterility	machine sterile
Floor outside	conical. pointed
Color graduation	blue
Format/ Packing	in bag
Material tube	Polypropylene (PP)
Material closure/ lid	Polypropylene (PP)
Division	0.5 ml
Type Tube/ Micro Vessel	Centrifuge tube
centrifugable up to	8500 g

