



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Nalgene centrifuge tubes PC, 28 ml, 25 x 92 mm,
10 pcs/pack**

**€1,042.00
plus VAT &
Shipping**

Product Images



Description

High speed polycarbonate centrifuge tubes with polypropylene screw cap.

- completely transparent
- round-bottomed
- loadable up to 50.000 x g in refrigerated and non-refrigerated centrifuges
- autoclavable
- temperature resistant from -135°C to +121°C

Additional Information

No.	CD-0763
Manufacturer (Brand)	Nalgene
old neolab article no.	2-9216
EAN	4058072219178
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	25 mm
Color	transparent (crystal clear)
Volume	28 ml
sterile	No
Height	92 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-135 °C
Floor outside	around
Format/ Packing	single. in carton/ dispenser
Material tube	Polycarbonate (PC)
Material closure/ lid	Polypropylene (PP)
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	50000 g

