



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene centrifuge tubes PC, 50 ml, 29 x 104 mm, 10 pcs/pack

**€151.10** As low  
as  
**€121.60**  
plus VAT &  
Shipping

### Product Images



## Description

---

High speed polycarbonate centrifuge tubes with polypropylene screw cap.

- completely transparent
- round-bottomed
- loadable up to 50.000 x g in refrigerated and non-refrigerated centrifuges
- autoclavable
- temperature resistant from -135°C to +121°C

## Additional Information

No.	CD-0765
Manufacturer (Brand)	Nalgene
old neoLab article no.	2-9217
EAN	4058072167905
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	29 mm
Color	transparent (crystal clear)
Volume	50 ml
sterile	No
Height	104 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-135 °C
Floor outside	around
Format/ Packing	single. in carton/ dispenser
Material tube	Polycarbonate (PC)
Material closure/ lid	Polypropylene (PP)
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	50000 g

