



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene centrifuge tubes PC, 85 ml, 38 x 106 mm, 10 pcs./pack

€171.50 As low
as
€138.10
plus VAT &
Shipping

Product Images



Description

High speed polycarbonate centrifuge tubes with polypropylene screw cap.

- completely transparent
- round-bottomed
- loadable up to 50.000 x g in refrigerated and non-refrigerated centrifuges
- autoclavable
- temperature resistant from -135°C to +121°C

Additional Information

No.	CD-0766
Manufacturer (Brand)	Nalgene
old neoLab article no.	2-9218
EAN	4058072219185
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	38 mm
Color	transparent (crystal clear)
Volume	85 ml
sterile	No
Height	106 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-135 °C
Floor outside	around
Format/ Packing	single. in carton/ dispenser
Material tube	Polycarbonate (PC)
Material closure/ lid	Polypropylene (PP)
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	50000 g

