



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene centrifuge tubes 15 ml, conical, made of PC, 10 pcs/pack

**€513.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Conical centrifuge tubes with exceptionally high mechanical strength.

- usable up to 6000 x g in refrigerated and non-refrigerated centrifuges
- co-molded graduation
- ideally suited for the separation of intact cells
- autoclavable
- crystal clear

## Additional Information

No.	CD-0737
Manufacturer (Brand)	Nalgene
old neolab article no.	2-9225
EAN	4058072219123
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	17 mm
Color	transparent (crystal clear)
Volume	15 ml
sterile	No
Height	100 mm
TBST MAX	121 °C
TBST MIN	-135 °C
Floor outside	conical
Color graduation	transparent (milky)
Format/ Packing	single. in carton/ dispenser
Material tube	Polycarbonate (PC)
Division	0.1 ml
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	6000 g

