



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Centrifuge tubes (PP) 13 ml, 16x100
mm, cap blue, 250 pcs./pack**

**€35.30
plus VAT &
Shipping**

Product Images



Description

Reagent and centrifuge tubes made of selected, biologically inert polypropylene.
With reliably sealing, blue PP screw cap, assembled.
For safe transport and storage of sample material.

Additional Information

No.	2-0140
Manufacturer (Brand)	neoLab
EAN	4058072047542
autoclavable at 121°C	No
diameter	16 mm
Color	transparent (milky)
Volume	13 ml
sterile	No
Height	100 mm
Type closure/ lid	Screw cap
TBST MAX	110 °C
TBST MIN	-10 °C
Floor outside	around
Format/ Packing	in bag
Material tube	Polypropylene (PP)
Material closure/ lid	Polypropylene (PP)
Type Tube/ Micro Vessel	Centrifuge tubes test tubes/tubes
Wall thickness	normal

