



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene centrifuge tubes PPCO 28 ml, 25x94 mm, screw cap PP, 10 pcs./pack

€106.60 As low
as
€85.80
plus VAT &
Shipping

Product Images



Description

High speed centrifuge tubes made of PPCO with screw cap made of PP.

Polypropylene copolymer can be autoclaved repeatedly and has better chemical resistance than PC.

- loadable up to 50.000 x g in refrigerated centrifuges
- round-bottomed
- leak-proof
- transparent

Temperature resistant: -40°C to +121°C

Additional Information

No.	CD-0768
Manufacturer (Brand)	Nalgene
old neolab article no.	2-9221
EAN	4058072219208
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	25 mm
Color	transparent (milky)
Volume	28 ml
sterile	No
Height	94 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-40 °C
Floor outside	around
Format/ Packing	single. in carton/ dispenser
Material tube	Polypropylene copolymer (PPCO)
Material closure/ lid	Polypropylene (PP)
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	50000 g

