



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene polysulfone centrifuge tubes, 30 ml, 26 x 95 mm, 10 pcs/pack

€1,265.00
plus VAT &
Shipping

Product Images



Description

Centrifuge tubes made of polysulfone with screw cap.

- transparent, slightly yellowish
- round bottom
- loadable up to 50.000 g
- temperature resistant from -100°C to +150°C
- resistant to acids and bases
- sterilizable with slightly attached lid

Additional Information

No.	CD-0753
Manufacturer (Brand)	Nalgene
old neoLab article no.	1-0031
EAN	4058072219147
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	25,3 mm
Color	transparent (orange)
Volume	30 ml
sterile	No
Height	92 mm
Type closure/ lid	Screw cap
TBST MAX	150 °C
TBST MIN	-100 °C
Floor outside	around
Format/ Packing	single. in carton/ dispenser
Material tube	Polysulfone (PSU)
Material closure/ lid	Polypropylene (PP)
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	50000 g

