



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoCulture Centrifuge Tubes PP 50 ml, with standing rim, (20 x 25) 500 pcs/pack, sterile

Special Price

€125.00 was

€162.40

**plus VAT &
Shipping**

Product Images



Description

Centrifuge tubes with white labeling field in excellent transparency and with high chemical resistance made of PP. With precise graduation by black imprint.

Additional Information

No.	C-8217
Manufacturer (Brand)	neoCulture
EAN	4058072154189
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	30 mm
Color	transparent (clear)
Volume	50 ml
sterile	Yes
Height	121 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-80 °C
Floor outside	conical. pointed . with stand ring
Color graduation	black
Type of sterility	radiation sterile
Division	2.5 ml
Material tube	Polypropylene (PP)
Material closure/ lid	High-density polyethylene (HDPE)
Type Tube/ Micro Vessel	Centrifuge tube
Wall thickness	normal
centrifugable up to	6000 g

