



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Ratiolab centrifuge tubes 50 ml, conical, sterilized

Special Price

€53.90 was

€82.00

**plus VAT &
Shipping**

Product Images



Description

Centrifuge tubes made of polypropylene, with conical bottom or with standing rim, in sterilized design. The sterile centrifuge tubes are certified RNase, DNase, endotoxin, pyrogen and ATP-free. The centrifuge tubes are provided with red screw cap, printed graduation and title block. Temperature resistance -80°C to 121°C

Additional Information

| | |
|-----------------------|---------------------|
| No. | RB-0040 |
| Manufacturer (Brand) | Ratiolab |
| EAN | 4058072984373 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | transparent (clear) |
| Volume | 50 ml |
| sterile | Yes |
| Type closure/ lid | Screw cap |
| TBST MAX | 121 °C |
| TBST MIN | -80 °C |
| Floor outside | conical |
| centrifugable up to | 7600 g |

