



neolab Migge GmbH

Rischerstr. 7-9

69123 Heidelberg

Deutschland

+49 (0)6221 / 8442-44

<https://www.neolab.de>

Umsatzsteuer-Identifikationsnummer:

DE 143 450 657



## Ratiolab centrifuge tubes 50 ml, standing rim, sterilized

Special Price

€53.20 was

€76.06

plus VAT & Shipping



## Description

Centrifuge tubes made of polypropylene, with conical bottom or with standing rim, in sterilized design. The sterile centrifuge tubes are certified RNase, DNase, endoxin, pyrogen and ATP-free. The centrifuge tubes are provided with red screw cap, printed graduation and title block. Temperature resistance -80°C to 121°C

## Additional Information

|                       |                          |
|-----------------------|--------------------------|
| No.                   | RB-0042                  |
| Manufacturer (Brand)  | Ratiolab                 |
| EAN                   | 4058072984397            |
| Transport temperature | Room temperature         |
| autoclavable at 121°C | Yes                      |
| Color                 | transparent (clear)      |
| Volume                | 50 ml                    |
| sterile               | Yes                      |
| Type closure/ lid     | Screw cap                |
| TBST MAX              | 121 °C                   |
| TBST MIN              | -80 °C                   |
| Floor outside         | conical. with stand ring |
| centrifugable up to   | 7600 g                   |

